

flexible

compact

mobile

# Emergency-ECG

**BIOSCOPE C**

## BIOSCOPE C - Pocket-Size ECG Monitor

Range of application: Emergency medicine - Out-patient visits - Pacemaker check-up



### BIOSCOPE C product features

Easy operation mode

- ECG detection on patients` chest by integrated 3-point-electrodes
- Cable adapter for connection of ECG patient cable
- ECG recording up to 11 minutes
- EKG view by internal LCD display
- Optional Analog- or Digital-adapter for ECG datatransfer to external ECG device or PC
- ECG-Analysis by Software BIOTREND

### Technical data BIOSCOPE C

- Lightweight, portable, battery powered (9V)
- Weight 350g - Size: 152 x 70 x 40 mm
- LCD display for ECG und Heartrate
- ECG recording up to 11 minutes
- 3 recording speeds (25,50, 100 mm/s)
- Autom. amplification adjustment (0,25...8 mV)
- Display of ST-Depression
- Detection and display of pacemaker impulses
- ECG-Patientcable for disposabel electrodes
- Option: Interface for connection to ECG device
- Option: Interface for connection to PC incl Software BIOTREND for ECG analysis
- Bandwith 0,5 ... 40 Hz (3db)
- Input impedance 2 M $\Omega$
- Defibrillator protection



### **Optional Accessories:**

- Cable adapter and ECG patient cable
- Analog-Interface for recording the stored ECGs
- Digital-Interface for connection PC incl. Software BIOTREND I.



MESA Medizintechnik GmbH  
Schärflmühlweg 4  
83671 Benediktbeuern  
Tel.: +49 (0) 8857 – 6918-0  
Fax: +49 (0) 8857 – 6918-29  
Email: info@mesamed.de  
www.mesamed.de

Our products are up from development to manufacturing subject to a certified Quality Management System in accordance to the EN ISO directives. The company is EN ISO certified. All products are labeled with CE mark and in accordance to the Medical Device Directives MDD 93/42/EEC.