

flexible

compact

# FOREHEAD-THERMOMETER

mobile

hygienic

**THERMOFINDER FS-300**

## THERMOFINDER FS-300 contactless forehead thermometer

Compact, small forehead thermometer for fast measurement of the body temperature without skin contact.



## THERMOFINDER FS-300

### Product features:

- Easy to operate, pocket size
- Fast non-contact measurement of body temperature by integrated infrared sensor (IR)
- Mode selection:; Body temperature and wide temperature (environment)
- LCD-Display: Temperature and battery capacity
- Data storage: Automatically up to 32 measurements, multi measurement
- Automatic power off function



### Technical data:

#### THERMOFINDER FS-300

- Light weight, compact (2 AAA batteries)
- Weight: 120 g - Size:62 x 42 x 103 mm
- LCD-Display (back lighted)
- Resolution and accuracy: 0,1 °C
- 5000 measurements per battery set
- Keys: Mode, start, storage, tone
- Measuring time: < 2 sec.
- Automatic power off (1 min.)
- Incl. Desk holder and batteries



MESA Medizintechnik GmbH  
Schärfilmühlweg 4  
D-83671 Benediktbeuern  
GERMANY  
Tel.: +49-8857 – 6918-0  
Fax: +49-8857 – 6918-29  
Email: [info@mesamed.de](mailto:info@mesamed.de)  
[www.mesamed.de](http://www.mesamed.de)

Our products are up from development to manufacturing subject to a certified Quality Management System in accordance to the EN ISO directives. The company is EN ISO certified. All products are labeled with CE mark and in accordance to the Medical Device Directives MDD 93/42/EEC.