

effective

SpO<sub>2</sub>

innovative

Oxygen saturation

# Pulse oximeter

light weight

compact

mobile

**Pulse oximeter – Finger pulse oximeter**

### Fingertip Pulse oximeter MS 20



- Digital display of heartrate and oxygen saturation
  - Comfortable application
  - Lightweight of only 50 g - robust housing
  - Automatic switch OFF function
  - operation possible up to 30 hours by use of 2 pcs. (AAA) batteries
  - Large LED display
- Easy operable Finger pulse oximeter with digital display of heartrate and SpO2  
Low cost version for fast and easy measurement

### Fingertip Pulse oximeter MS 300



- Graphic display of heartrate as diagram
  - 2 colour display with 10 programmable levels of lightness
  - 6 programmable display views
  - Lightweight of only 50 g - robust housing
  - Automatic switch OFF function
  - Betrieb bis zu 30 Stunden mit 2 x herkömmlichen (AAA) Batterien
- The digital Finger pulse oximeter MS 300 is an oxygen and heartrate measuring device for hospital, rescue service and home care  
Programmable lightness and different display views  
Documentation as bar or wave diagram for optimal operation conditions  
Small and light weight and designed for operation during night time

### Pulse oximeter PM50



- For all patient groups (adult, pediatric, neonatal)
- For all clinical ranges f.e. rescue services
- Large numeric display of characters
- High resolution LCD display for trend and Spo2 data
- Automatic Standby and Power-off function
- Up to 100 patienten-ID's and 200 measurements
- Nellcor Oxisensor compatible

### Pulse oximeter PM60



- For all patient groups (adult, pediatric, neonatal)
- For all clinical ranges f.e. rescue services
- Large numeric display of characters
- High resolution colour LCD display for trend and Spo2 data
- Acoustic and visual alarm
- Pitch tone
- Infrared technology
- Programmable alarm tones / alarm tone limits
- 96 h trend memory up to 4000 Spotcheck measurements
- Nellcor Oxisensor compatible

Scope of delivery MS20, MS300: Fingertip pulse oximeter, 2 AAA batteries, clip band, manual  
Scope of delivery PM50: Pulse oximeter PM50, 4 AA batteries, Fingerclip sensor adult, holding band, manual  
Scope of delivery PM60: Pulse oximeter PM60, 3 AA batteries, Fingerclip sensor adult., holding band, manual



MESA Medizintechnik GmbH  
Schärflmühlweg 4  
D-83671 Benediktbeuern  
GERMANY  
Tel.: +49 (0) 8857 – 6918-0  
Fax: +49 (0) 8857 – 6918-29  
Email: info@mesamed.de  
www.mesamed.de

Our products are up from development to manufacturing subject to a certified Quality Management System in accordance to the EN ISO directives. The company is EN ISO certified. All products are labeled with CE mark and in accordance to the Medical Device Directives MDD 93/42/EEC.