

compact
qualitative
functional
Electrode suction system
stable
flexible
mobile

DT 80 B Series
DT 80 T Series

DT 80 B and DT 80 T vacuum electrode suction systems for effective ECG recordings are available as different models.

Electrode suction devices DT 80 B and DT80 T :

DT 80 B and the DT 80 T vacuum electrode suction systems perform high quality ECG recordings in a fast, effective and also comfortable way. These specific electrodes suction units are compatible with nearly all ECG devices and can be installed onto various trolley types, too.

The DT 80 B solid turnaround holder arm and the DT 80 T telescope holder arm are constructed for heavy duty and sophisticated use. Quality and performance are outstanding characteristics of these systems. An adjustable brake for the vertical flexibility is integrated additionally.

DT 80 B & DT 80 T features:

- Flexible holder arm with integrated defibrillator protection
- DT 80 B arm turnaround range 315°, 120/240 cm
- DT 80 T Telescope arm, turnaround range 360°, vertical radius 150°, length variable from 84 cm up to 132 cm
- Easily to operate
- Vacuum pump with blow out function and four level adjusting
- Automatic regulation of vacuum power for every level
- Vacuum range 60 – 200 mbar, display of pressure status
- Switchable pump function for end of suction and with blow out function for cleaning of electrodes
- High safety by protected currency connector
- High qualitative ECG recordings by shielded leads
- Single changeable suction leads (plug in principle)
- Ag/AgCl suction electrodes guarantee optimal adhesive force, even at patients with strong hair growth
- Different lengths of suction leads (10x1m or 6x1m & 4x1,3m)



DT 80 B plus

DT 80 T plus

The DT 80B and DT 80 T electrode suction systems can be mounted to different trolleys. Additionally desk- and wallholder models are available. The MESA MW-1250-60 trolley and other trolleys are designed for stationary or mobile use. There are available as several models and flexible concerning dimensions of e.g height, bridth, and number of plates.



DT 80 B system DT 80 TH system
both with trolley H



DT 80 B system DT 80 T system
both for table mounting



DT 80 T MW system
Mounting model with MW-1250-60 trolley



MESA Medizintechnik GmbH

Schärfelmühlweg 4

D-83671 Benediktbeuern

GERMANY

Tel.: +49 (0) 8857 – 6918-0

Fax: +49 (0) 8857 – 6918-29

Email: info@mesamed.de

www.mesamed.de

Our products are up from development to manufacturing subject to Quality Management System. The company is EN ISO certified. All products are labelled with CE mark.

The illustrations within this brochure show some optional features. All information's represent the technical product status at date of printing. Changes due to technical advancement are subject to modification of amendment. Manufacturer: Strässle & Co. Medizintechnik GmbH (Rev. A03-2023).