

flexible

compact

mobile

Emergency-ECG

ECG Monitor PC-80 B Easy ECG

PC-80 B Easy ECG – Pocket size ECG Monitor with LCD colour display

Compact and small ECG-Monitor with memory function for fast recording of an emergency ECG



PC-80 B Easy ECG produkt features

- Easy operation mode
- Fast ECG recording by contact with integrated ECG-Elektrodes at Finger and palm of hand or at finger and chest
- Connection for 3-lead ECG patient cable
- ECG recoding and memory up to 1200 ECGs
- Up to 10 h ECG recording with one battery set
- Indication of capacity of batteries
- ECG view by internal LCD colour display
- USB cable for data transfer to PC
- Software for analysis and documentation
- Delivery with ECG cable, pouch, manual, Software on CD and USB cable

Technical data PC-80 B Easy ECG

- Light, compact, battery powered (2 AAA)
- Weight 110 g, Size: 125 x 70 x 21,5 mm
- LCD colour display for ECG and Heartrate (57,6 x 40 mm)
- ECG detection rate: 1 Hz – 40 Hz
- ECG alarm tone
- Heartrate: 30 bpm up to 240 bpm
- Recording (Monitoring) / ECG max.30 sec.
- Memory capacity up to 200 ECG recordings
- ECG analysis within 30 sec.
- PC connection incl. documentation Software
- EMC Group I, Class B
- Protection class B



MESA Medizintechnik GmbH
Schärflmühlweg 4
83671 Benediktbeuern
Tel.: +49 (0) 8857 – 6918-0
Fax: +49 (0) 8857 – 6918-29
Email: info@mesamed.de
www.mesamed.de

Our products are up from development to manufacturing subject to a certified Quality Management System in accordance to the EN ISO directives. The company is EN ISO certified. All products are labeled with CE mark and in accordance to the Medical Device Directives MDD 93/42/EEC.