

effective

sturdy

# Centrifuge

low noise level

compact

**SERVOspin III**



Equipped with a 8 - angle rotor for tubes with a volume of 15 ml and a max. RCF of 3.074 g the SERVOspin III is the optimal device for a general practitioner. The SERVOspin III centrifuge needs just little storage place and is very easy to use. The last selected setting for rotational speed and time can be saved. It is ready to use at the touch of the button. For users safety, the lid gets locked electromagnetically during operation. The durable plastic case and hygienic membrane keyboard enable an easy cleaning.



**Technical Data SERVOspinIII:**

Max. speed : 100 - 5000 r/min

Rotation frequency: increase by 100 r/min

Max.RCF (relative centrifugal force): 3074 G

Rotor: 8 x 15 ml

Noise: app. 56 dB (A)

Unbalance control: yes

Rotations / RCF: yes

Dimensions: 350 x 400 x 320 mm

Weight: 5 kg

Power: 220 - 230 Volt 50/60 Hz

Inner  $\varnothing$  of insets x length: 15 x 105 mm



MESA Medizintechnik GmbH  
Schärflmühlweg 4

D-83671 Benediktbeuern  
GERMANY

Tel.: +49 (0) 8857 – 6918-0

Fax: +49 (0) 8857 – 6918-29

Email: [info@mesamed.de](mailto:info@mesamed.de)

[www.mesamed.de](http://www.mesamed.de)

Our products are up from development to manufacturing subject to a certified Quality Management System in accordance to the EN ISO directives. The company is EN ISO certified. All products are labeled with CE mark and in accordance to the Medical Device Directives MDD 93/42/EEC.