

effective

sturdy

# Centrifuge

low noise level

compact

**SERVOspin III**



Equipped with a 8 - angle rotor for tubes with a volume of 15 ml and a max. RCF of 3.074 g the SERVOspin Plus is the optimal device for a general practitioner, a hospital or a laboratory. The SERVOspin Plus centrifuge needs just little storage place and is very easy to use. The last selected setting for rotational speed and time can be saved. It is ready to use at the touch of the button. For the safety of the user the lid gets locked electromagnetically during operation. The durable plastic case and hygienic membrane keyboard enable an easy cleaning procedure.



#### **Technical Data SERVOspin Plus:**

Max. speed : 100 - 5000 r/min  
Rotation frequency: increase by 100 r/min  
Max.RCF (relative centrifugal force): 3074 G  
Rotor: 8 x 15 ml  
Noise: app. 56 dB (A)  
Unbalance control: yes  
Rotations / RCF: yes  
Dimensions: 350 x 400 x 320 mm  
Weight: 5 kg  
Power: 220 - 230 Volt 50/60 Hz  
Inner  $\varnothing$  of insets x length: 15 x 105 mm



MESA Medizintechnik GmbH  
Schärflmühlweg 4  
D-83671 Benediktbeuern  
GERMANY  
Tel.: +49 (0) 8857 – 6918-0  
Fax: +49 (0) 8857 – 6918-29  
Email: [info@mesamed.de](mailto:info@mesamed.de)  
[www.mesamed.de](http://www.mesamed.de)

Our products are up from development to manufacturing subject to a certified Quality Management System in accordance to the EN ISO 13485:2016 directives. All products are CE marked.