



[www.mesamed.de](http://www.mesamed.de)

## DRAEGER MULTIMED POD MULTIPARAMETER SYSTEM ECG - SPO2 - TEMPERATURE - NIBP

### Multiparameter - Trunk cable



#### **P/N MON-801530**

Multimed Pod Monitoring Multiparameter-Trunk cable  
for SC 6000/7000/9000 - Draeger Ref. 3368391  
(compatible to Original Draeger cable holder)



#### **P/N ZUB-801530**

Holder for the multi-parameter cable yoke  
P/N MON-801530 and Draeger Ref. 3368391  
(to be fixed onto the pole or wall on bedside)



[www.mesamed.de](http://www.mesamed.de)

## DRAEGER MULTIMED POD MULTIPARAMETER SYSTEM ECG - SPO2 - TEMPERATURE - NIBP

### ESU-Block



**P/N DRÄ-5947226**

Draeger ECG ESU-Block 5-lead (HF-Filter)

Draeger Ref. 5947226

to be connected to ECG leadwire sets

P/N **ADK-101430** or

P/N **ADK- 101450** or

**The ESU block is not required for the following ECG cable:**

P/N **ADK-101656** and

P/N **ADK-101653**



[www.mesamed.de](http://www.mesamed.de)

## DRAEGER MULTIMED POD MULTIPARAMETER SYSTEM ECG - SPO2 - TEMPERATURE - NIBP

### ECG patient leadwire sets for Multimed Pod



#### **P/N ADK-101430**

Set of 3-lead ECG leadwires (onto 5-lead block),  
with clip, VS block system, length 100 cm, shielded,  
colour code: red, green, yellow

for Multimed Pod – Draeger Ref. MP03411

OR-Version with ESU-Filter **P/N ADK-101656** (length 90 cm)



#### **P/N ADK-101440**

Set of 5-lead ECG leadwires with clip, VS block system,  
length 100 cm, shielded,

colour code: black, white, red, green, yellow  
for Multimed Pod – Draeger Ref. MP03413


OR-Version with ESU-Filter **P/N ADK-101653** (length 90 cm)



[www.mesamed.de](http://www.mesamed.de)

## DRAEGER MULTIMED POD MULTIPARAMETER SYSTEM ECG - SPO2 - TEMPERATURE - NIBP

### SpO2-Adapter cable Multimed Pod



#### **P/N SPA-7512**

SpO2 adapter cable Multimed Pod  
for Monitors SC6000/7000/9000,  
length 1 m,  
Draeger Ref. 3368433



#### **P/N SPA-7513**

SpO2 adapter cable Multimed Pod  
for Monitors SC6000/7000/9000  
length 2 m,  
Draeger Ref. 3375834



[www.mesamed.de](http://www.mesamed.de)

## DRAEGER MULTIMED POD MULTIPARAMETER SYSTEM ECG - SPO2 - TEMPERATURE - NIBP

### SpO2-Sensors



**P/N FM-6212-90**

Fingerclip sensor Type FM-6212-90  
Nellcor Oxisensor compatible,  
length 90 cm,  
to be connected to SpO2 adapter cable  
**P/N SPA-7512 or SPA-7513**



**P/N SPR-3212-9 plus**

SpO2 SoftTip Plus Sensor Type R-3212-9 plus,  
large, for adult,  
Nellcor Oxisensor compatible (autoclavable)  
Length: 90 cm  
to be connected to SpO2 adapter cable  
**P/N SPA-7512 or SPA-7513**

**Other Original Nellcor Oximax und Nellcor Oxisensor compatible SpO2-Sensors can be found at**

[http://www.mesamed.de/mesa-ftp/pdf/ger/Katalog\\_Nellcor.pdf](http://www.mesamed.de/mesa-ftp/pdf/ger/Katalog_Nellcor.pdf)


Nellcor is a registered trademark of Nellcor Puritan Bennett Inc..



[www.mesamed.de](http://www.mesamed.de)


## DRAEGER MULTIMED POD MULTIPARAMETER SYSTEM ECG - SPO2 - TEMPERATURE - NIBP

### SpO2-disposable-Sensors



#### **P/N SPD-2211-1**

SpO2-disposable-Sensors Typee DA-2211-1 Nellcor Oxisensor,  
Length: 45 cm / Box of 24 peace, Micro-Foam,  
adult - Patientweight: >30 kg,  
to be connected to SpO2-Adapter cable **P/N SPA-7512 or SPA-7513**



#### **P/N SPD-2211-2**

SpO2-disposable-Sensors Typee DA-2211-2 Nellcor Oxisensor  
Length: 45 cm / Box of 24 pcs., Micro-Foam,  
Pediatric - Patientweight: 10 – 50 kg,  
to be connected to SpO2-Adapter cable **P/N SPA-7512 or SPA-7513**

**More compatible Spo2 sensors e.g.  
SoftTip Plus medium, SoftTip Plus small,  
MultiSoft, FlexMon, Disposable SpO2-Sensors,  
Earsensors in various versions can be found there**

[http://www.mesamed.de/mesa-ftp/pdf/ger/Katalog\\_Sensoren.pdf](http://www.mesamed.de/mesa-ftp/pdf/ger/Katalog_Sensoren.pdf)

**There you can also find many more Spo2-AdapterCable,  
SpO2-Sensors, Oxygensensors and Flowsensors  
compatible to many manufacturers**



[www.mesamed.de](http://www.mesamed.de)

## DRAEGER MULTIMED POD MULTIPARAMETER SYSTEM ECG - SPO2 - TEMPERATURE - NIBP


### Temperature Probes

Connection cable for disposable Temperature Probes



#### **P/N TMP-120010**

Connection cable for disposable temperature probes  
Type: DRAEGER/Siemens, Length 3m



#### **P/N TMP-120011**

Compatible disposable temperature  
probe, skin contact, length 75 cm  
1 Box = 20 pcs,

#### **P/N TMP-120012**

Compatible disposable temperature  
probe, rectal-oesophageal, length 90 cm,  
1 Box = 20 pcs,


To be connected to adapter cable **P/N TMP-120010**



[www.mesamed.de](http://www.mesamed.de)

## **DRAEGER compatible ECG-Trunkcable and ECG-leads, Spo2- Adaptercable and Spo2-Sensors**

Reusable temperature probes with direct connect



### **P/N SI-SAR-275A**

Temperature Probe adult, skin,  
reusable, Length: 1,5 m, Ø 18 mm, autoclavable,  
Type SIEMENS Ref. 4329822



### **P/N SI-PR-275A**

Temperature probe pediatric, skin,  
reusable, Length: 2,75 m, Ø 10 mm, autoclavable,  
Type SIEMENS



### **P/N SI-NR-275A**

Temperature probe Infant/Neonatal, skin,  
reusable, Length: 2,75 m, Ø 6 mm, autoclavable,  
Type SIEMENS



**YSI compatible connection cables and reusable YSI compatible  
Temperature probes on request!**





[www.mesamed.de](http://www.mesamed.de)

## **DRAEGER compatible ECG-Trunkcable and ECG-leads, Spo2-Adaptercable and SpO2-Sensors**

Disposable-cover for reusable skin temperature probes



### **P/N TMP-150010**

Disposable temperature probe cover  
for use with reusable skin temperature  
probes to be fixed onto the skin and  
reflect the ambient light to ensure accurate  
measurement

1 Box = 12 pcs.



[www.mesamed.de](http://www.mesamed.de)

## **DRAEGER compatible ECG-Trunkcable and ECG-leads, Spo2-Adaptercable and Spo2-Sensors**

Reusable temperature probes with direct connect



### **P/N SI-AR-275A**

temperature probe adult, rectal,  
reusable, Length: 2,75 m, Ø 4,8 mm, autoclavable,  
Type SIEMENS - Ref. 4329889



### **P/N SI-PR-275A**

temperature probe pediatric, rectal,  
reusable, Length: 2,75 m, Ø 4,0 mm, autoclavable,  
Type SIEMENS



### **P/N SI-NR-275A**

temperature probe Infant/Neonatal, rectal,  
reusable, Length: 2,75 m, Ø 2,7 mm, autoclavable,  
Type SIEMENS



**YSI compatible connection cables and reusable YSI compatible temperature probes on request!**



[www.mesamed.de](http://www.mesamed.de)

## DRAEGER MULTIMED POD MULTIPARAMETER SYSTEM ECG - SPO2 - TEMPERATURE - NIBP

### NIBP Air hose



Monitor side

Cuff side

#### **P/N MAN-130010**

NIBP-Air hose, length 2,5 m, TPU,  
Single tube system for adult and pediatric  
Type: Philips A1/A3, Nellcor, DRAEGER/Siemens,  
Datascope, Spacelabs, Kontron, Mindray



[www.mesamed.de](http://www.mesamed.de)

## DRAEGER MULTIMED POD MULTIPARAMETER SYSTEM ECG - SPO2 - TEMPERATURE - NIBP

Disposable NIBP cuffs



### **P/N MAN-130015**

NIBP-Disposable-Cuff, non-woven,  
white, 27,5 - 36,5 cm,

Single tube system for adult

Type: Philips, DRAEGER/ Siemens, Datascope, Spacelabs

1 Box = 5 pcs.

to connect with air hose P/N **MAN-130010**

### **P/N MAN-130016**

NIBP-Cuff, non-woven,  
white, 13,8 - 21,5 cm,

Single tube system for pediatric,

Type: Philips, DRAEGER/ Siemens, Datascope, Spacelabs

1 Box = 5 pcs.

to connect with air hose P/N **MAN-130010**



[www.mesamed.de](http://www.mesamed.de)

## DRAEGER MULTIMED POD MULTIPARAMETER SYSTEM ECG - SPO2 - TEMPERATURE - NIBP

NIBP cuffs, reusable, TPU



### **P/N MAN-130013**

NIBP-Cuff, reusable, TPU,  
25 - 35 cm, grey

Single tube system for adult

Type: Philips (HP), Siemens/DRAEGER, Datascope, Spacelabs,  
to connect with air hose P/N **MAN-130010**

### **P/N MAN-130014**

NIBP-Cuff, reusable, TPU,  
33 - 47 cm, violet

Single tube system for adult (large)

Type: Philips (HP), Siemens/DRAEGER, Datascope, Spacelabs,  
to connect with air hose P/N **MAN-130010**



[www.mesamed.de](http://www.mesamed.de)

## DRAEGER MULTIMED POD MULTIPARAMETER SYSTEM ECG - SPO2 - TEMPERATURE - NIBP

NIBP-Cuffs, reusable, TPU,  
Colour codes:

Adult	
Adult long	
Small Adult	
Pediatric	
Infant	
Neonate	
Large Adult	
Large adult long	
Thigh	

Rev.A01-2012



TELEFON (08857) 6918- 0

TELEFAX (08857) 6918-29

Schärfmühlweg 4

**D-83671 Benediktbeuern**

[info@mesamed.de](mailto:info@mesamed.de)