

compact

qualitative


functional

# Enviromental Parameter Analyzer

weather station

**EPA-3**





## SPECIFICATION

### Measurement range

Temperature	0-50 °C
Relative humidity	0-100%
Air pressure	750-1200 hPa

### Accuracy

Temperature	±1%
Relative humidity	±3%
Air pressure	±1%

<b>Communication</b>	USB
----------------------	-----



The EPA-3 Enviromental Parameter Analyzer is a full-featured weather station compatible with MESA PRO SPIRO spirometer products. EPA-3 can be connected to computers via USB and automatically will set the environmental conditions in the PRO SPIRO spirometer software..



MESA Medizintechnik GmbH  
Schärflmühlweg 4  
D-83671 Benediktbeuern  
GERMANY  
Tel.: +49 (0) 8857 – 6918-0  
Fax: +49 (0) 8857 – 6918-29  
Email: [info@mesamed.de](mailto:info@mesamed.de)  
[www.mesamed.de](http://www.mesamed.de)

Our products are up from development to manufacturing subject to a certified Quality Management System in accordance to the EN ISO 13485:2016 directives. All products are CE marked.