

effective

SpO₂

innovative

Oxygen saturation

Pulse Oximeter

light weight

compact

mobile

Pulse Oximeter

MESA provides a comprehensive Pulse oximeter program in different versions.



MS 20 Fingertip Pulse Oximeter, Rubberpouch, textile holster and clip band



MS 300 Fingertip Pulse Oximeter with integrated pulswave and programmable alarm tones, rubber pouch, textile holster and clip band



MS 800 Handheld Pulse Oximeter with integrated 3-Channel ECG programmable alarm tones, Li-Ion battery, Charging unit



MS 600 Handheld Pulse Oximeter with backlighted TFT-Colour display, integrated pulse wave, programmable alarm limits, profiles and alarm tones, various SpO2-Sensors available, optional: Temperature probe, rechargeable Li-Battery, PC-Software



MS 700 Handheld Pulse Oximeter with Touchscreen, integrated pulse wave, and programmable alarm limits and alarm tones, rechargeable LiPo battery, optional: PC-Software



MESA Medizintechnik GmbH
 Schärflmühlweg 4
 D-83671 Benediktbeuern
 GERMANY
 Tel.: +49 (0) 8857 – 6918-0
 Fax: +49 (0) 8857 – 6918-29
 Email: info@mesamed.de
www.mesamed.de

Our products are up from development to manufacturing subject to a certified Quality Management System in accordance to the EN ISO directives. The company is EN ISO certified. All products are labeled with CE mark and in accordance to the Medical Device Directives MDD 93/42/EEC.