

OxiPen Pulsoximeter

...inclusive Fingerclipsensor



Monitoring parameter:
Oxygen saturation, Heartrate

➤ **Included in delivery:**

- Portable Pulsoximeter OxiPen
- Fingerclipsensor SPF-3227
- 2x Battery Type AAA
- English manual

➤ **OPTION:**

- Y-Sensor Type SPY-3227
- Earsensor Type SPES-3227
- Small Fingerclipsensor Type SPFS-3227
- Wrapsensor Type SPW-3227
- Einmalsensoren Typ SPD-3227-*

Y-Sensor Type SPY-3226



Earsensor



Technical data:

- **SpO2**
 - **Measurement range:** 45-100%
 - **Accuracy:** +/- 2% at 70-100%
- **Pulse frequency**
 - **Meas.range:** 20 - 300 bpm
 - **Accuracy:** +/- 1 bpm up to 100 bpm and +/- 1% >100 bpm
- **LCD Display**
 - LCD backlighted with red, orange and green LED. View of %SpO2, Pulserate, Signalquality, Bargraph, Batterystatus and Sensorstatus
- **Pulseton**
 - Tone variable (depending on the measured O2- Saturation)
- **Signal**
 - Acoustic and visual Signals (Pulsefrequency, Sensorstatus, System failure and Batterystatus) Signal: ON/OFF function with automatic ON function after 3 minutes.
- **Power supply**
 - 2 x AAA ALKALINE
- **Operation time 15 hours**

Distributed by:



MESA Medizintechnik
GmbH
Postfach 1152
Schärflmühlweg 4
D-83671 Benediktbeuern
Tel.: +49 8857 6918-0
Fax: +49 8857 6918-29
www.mesamed.de
info@mesamed.de