

flexible

compact

INFRARED-THERMOMETER

mobile

hygienic

ThermoCheck TC 700 Forehead Thermometer



ThermoCheck TC 700 non-contact Forehead Infrared Thermometer

Compact and ergonomic forehead thermometer for fast measurement of the body temperature without skin contact. Body temperature measurement with fast results is possible for babies, children and adults. TC 700 avoids any risk of infections and guarantees operation by highest comfort.



TC 700 Non-Contact Forehead Thermometer - Product features:

- Easy to operate, pocket size
- Fast non-contact measurement of body temperature by infrared sensor (IR)
- Measuring distance 1 – 5 cm
- Mode selection:; Body temperature and temperature of environment
- LCD-Display: Temperature & battery capacity
- Data storage: Automatically up to 20 measurements, multi measurement
- Automatic power off function

Technical data ThermoCheck TC 700 Infrared Thermometer:

- Light weight and compact 150x88x40 mm
- Battery operated: (2 x AAA)
- HD-LCD-Display (back lighted)
- Resolution: +/- 0,1 °C, 0,1° F
- Accuracy: +/- 0,2 °C, 0,2° F
- Measuring range: 34,9° C up to 42,2 ° C
- Keys: Mode, start, storage, tone
- Measuring time: app. 1 sec.
- Automatic power off (after 10 +/- 1 sec.)
- Weight: 109,5 g (incl batteries)
- Incl. manual, 2 x AAAA batteries and textile bag



MESA Medizintechnik GmbH
Schärflmühlweg 4
D-83671 Benediktbeuern
GERMANY
Tel.: +49-8857 – 6918-0
Fax: +49-8857 – 6918-29
Email: info@mesamed.de
www.mesamed.de

Our products are up from development to manufacturing subject to a certified Quality Management System in accordance to the EN ISO 13485:2016 directives. All products are CE marked.

