

exact

Spirometer

Calibration syringue pump

mobile

compact

PCS 3000

PCS 3000 - the precise spirometer calibration syringe pump

The calibration syringe pump PCS 3000 can be used to check, adjust and calibrate all kinds of spirometers. The measuring volume is 3000 ml. This calibration pump is particularly suitable for on-site applications. The calibration pump is designed in a way that all spirometers can be directly connected for calibration. If the spirometer to be tested is connected via the flexible rubber connection piece, an immediate calibration is possible. The test volume of 3000 ml is generated with the hand lever by continuous pumping movements. The readings as reference for the function of the spirometer can be used to check the accuracy or to calibrate the spirometer.



MESA Medizintechnik GmbH
Schärflmühlweg 4
D-83671 Benediktbeuern
GERMANY
Tel.: +49 (0) 8857 – 6918-0
Fax: +49 (0) 8857 – 6918-29
Email: info@mesamed.de
www.mesamed.de

Our products are up from development to manufacturing subject to a certified Quality Management System in accordance to the EN ISO directives. The company is EN ISO certified. All products are labeled with CE mark and in accordance to the Medical Device Directives MDD 93/42/EEC.