



Capnograph CP 500



Capnograph CP 500 Type MEMO ETCO Mainstream

Small and super lightweight unit with ultimate portability.
Rapid monitoring of end tidal carbon dioxide (EtCO₂)
and respiration rate (RR) of adult, pediatric and neonatal patients.
Set includes 2 x airway adapter for adult/pediatric
and 2 x AAA battery.

Accessories:

airway adapter adult/pediatric
(disposable)



airway adapter neonatal
(disposable)



MESA Medizintechnik GmbH
Schärfmühlweg 4
D-83671 Benediktbeuern
Tel.: 0049-(0)8857 – 6918-0
Fax: 0049-(0)8857 – 6918-29
Email: info@mesamed.de
www.mesamed.de

