

compact
qualitative
functional
Medical Trolleys
stabile
flexible
mobile

Trolley-Program

Trolley MGW-series: MGW-1250-60 and MGW-1000-44

- High quality powder-coated working plates; aluminum profiles for assembly of ECG cable holder or suction devices
- Height adjustable working plates for Medical Devices e.g. ECG, PC or printer (color of plates RAL 7035)

Trolley MGW-1250-60

- Stationary and mobile trolley with aluminum profile frame and 4 plates, 2 plates with grip
- For assembly of ECG cable holder (optional) or electrode suction system by special holder for arm of suction device (optional)
- Extendable Mouse plate mounted on right side of 2nd upper working plate
- Mobile by 4 rolls, front rolls fixable, high stability
- Vacuum pump for electrode suction system to be placed on basement plate
- Dimensions: Trolley MGW-1250-60
125 x 60 x 44 cm (H x W x D),
weight: app. 14 kg



Trolley MGW-1000-44

- Stationary and mobile trolley with aluminum profile and grip on back side, 2 working plates
- For assembly of ECG cable holder (optional) or electrode suction system by special holder for arm of suction device (optional)
- Accessories basket flexible in height assembly
- Mobile by 4 rolls, front rolls fixable, good stability
- Vacuum pump for electrode suction system to be placed on basement plate
- Dimensions: Trolley MGW-1000-44
125 x 40 x 46 cm (H x W x D),
weight: app. 10 kg



MESA Medizintechnik GmbH
Schärflmühlweg 4
D-83671 Benediktbeuern
GERMANY
Tel.: +49 (0) 8857 – 6918-0
Fax: +49 (0) 8857 – 6918-29
Email: info@mesamed.de
www.mesamed.de