

Pulswave

SpO₂

Colour display

Oxygen saturation

Pulse oximeter+Capnograph

CO₂

mobile

compact

**Hand Held-Pulsoximeter MS 900
with integrated Capnograph**



The Hand Held Pulse Oximeter MS900 is a combination device as pulse oximeter with integrated capnograph. This device provides fast and reliable Heart rate (HR), Oxygen saturation (SpO2) and all relevant parameters of the Capnography (CO2). The bright 7.5 cm OLED color display shows the determined values easily to be detected.



Specifications of MS900 Type PC-900B Hand Held Pulse oximeter with Capnograph

- Clear 7.5 cm TFT color display with controllable brightness
- Temperature, pressure compensation and zero adjustment system
- Balanced gas compensation minimizes impact on other gases
- Lung pressure and temperature compensation (BTPS compensation) available
- Suitable for intubated or non-intubated patients
- Pump flow rate adjustable for sample lines for adults, small children and low flow samples
- Water filter especially for high humidity environments
- Adjustable audible and visual alarms for exceeding limit values and apnoea
- The 24-hour graphic trend can be saved after switching off the device
- MoveOxy SpO2 Fingerclip Sensor for Adult
- Integrated rechargeable Lithium-Ion batteries with up to 10 hours working time
- USB connectivity for data transfer
- Weight: 0.6 kg - Dimensions: 70 x 160 x 40 mm (WxHxD)



MESA Medizintechnik GmbH
Schärflmühlweg 4
D-83671 Benediktbeuern
Tel.: +49-8857 – 6918-0
Fax: +49-8857 – 6918-29
Email: info@mesamed.de
www.mesamed.de

Our products are up from development to manufacturing subject to a certified Quality Management System in accordance to the EN ISO directives. The company is EN ISO certified. All products are labeled with CE mark and in accordance to the Medical Device Directives MDD 93/42/EEC.